


Calibration Protocol with 3D Slicer:

Required Materials:


- Aurora EM Generator
 - Aurora EM Coils (x2) (inside probes)
 - Target Plate
 - Pen Probe
 - EIT Probe
 - Target Plate, Pen Probe, and EIT Probe STL files
 - 3D Slicer v5
 - Plus Server Launcher
-

Procedure:


1. Connect probes to Aurora and Aurora to the computer.
 - a. Connect Pen Probe to Port 1
 - b. Connect Probe Coil to Port 2
2. Turn on Aurora
3. Open 3D Slicer 5
4. Launch Plus Server Launcher
5. Set Device Configuration Directory to 'Aurora Two Tools'
6. Click Launch and should hear beep
7. Upload aurora_target_flex.STL, flexProbe.STL, and NeedleModel.STL files to 3D slicer
8. From 3D Slicer, go to **IGT** → **OpenIGTLinkIF**




9. Add Connector IO




10. Go to IGT → Pivot Calibration




11. Set Input and Output to *PenCoilToEM*



12. Place the Pen Coil on a point on the calibration plate
13. Click **Start Pivot Calibration**
14. Pivot the pen steadily multiple times.
15. Redo until $RMS \leq 0.1$
16. Go to **IGT** → **Fiducial Registration Wizard**




17. Mark the STL File at each of the four points
18. Place the pen in each of the four points marked
19. RMS <= .5 is acceptable
20. Place the probe in the calibration plate
21. Output the transformation matrix for SecondCoilToEM and CalPlate To EM



22. Create a new transformation called PCBtoSecondCoil
23. In MATLAB, calculate the PCBtoSecondCoil matrix from:

$T_{\text{plate_to_coil}} = \text{inv}(T_{\text{coil_to_em}}) * T_{\text{plate_to_em}};$

24. Copy the MATLAB output and edit PCBtoSecondCoil transformation to be the MATLAB output



25. The final transformation hierarchy should be as such:

